


### Closed Suction Catheter – Double Swivel Elbow



#### Description :

1	Cap for suction side
2	Control button
3	Sealed Sleeve
4	Suction Tubing
5	Atraumatic Suction Hole
6	Ventilator Port
7	Patient Port
8	Irrigation Port
9	MDI Port

#### For Endotracheal

Length	54 cm (21.3 in.)
Diameter	6Fr to 16 Fr
MDI	W/O

#### For Tracheostomy

Length	30.5 cm (12 in.)
Diameter	6Fr to 16 Fr
MDI	W/O

#### Instruction For Use:

1. Select appropriate size closed suction catheter.
2. Attach thumb control valve to suction tubing.
3. Depress and hold thumb valve and simultaneously adjust vacuum regulator to desired level.
4. Release thumb control valve and attach catheter between patient and the ventilator circuit.

#### Suggested suction procedure:

1. Stabilize catheter and Endotracheal adapter with one hand then push the catheter into the Endotracheal tube with the catheter into the Endotracheal tube with the thumb and forefinger of the opposite hand.
2. Advance catheter to desired depth.
3. Depress and hold thumb control valve, then gently withdraw catheter, **stop withdrawal when black marking ring is visible inside sleeve.**
4. Release thumb control valve
5. Repeat steps 1-4 above as necessary.