


# UROLOGY & GYNAECOLOGY



Adress: Loujin Industrial park, Shuofang, Wuxi, 214143, P.R.China



Tel: 86 510 8273 3513 Fax: 86 510 8273 0717



E-mail: [info@gcmedica.com](mailto:info@gcmedica.com)



Website: [www.gcmedica.com](http://www.gcmedica.com)



GCMEDICA ENTERPRISE LTD.,(WUXI )



## GCMEDICA PROFILE

GCMEDICA Group founded in 1998 with a mission dedicating to all medical personnel and patients, combines many years of experience and expertise, GCMEDICA specializes in the manufacturing scientific research and development of disposable medical devices and advanced moulds. GCMEDICA have built up a leading group company concentrated in various products production with CE certificate and registered with FDA US. With persistent and diligent exploration, GCMEDICA extents the product ranges to Respiratory, Anesthesia, Urology, Surgery, Operating room, Gynecology products.

GCMEDICA located in new industrial zone-Wuxi, which owns excellent convenient transportation. GCMEDICA has building area of 10000m<sup>2</sup> manufacture base which incorporates a 100,000-class clean room of 5500m<sup>2</sup> and well-equipped laboratories. GCMEDICA mainly manufactures Class I, Class II and Class III disposable medical devices. The applied expertise of our R&D team ensures our latest products are always innovative.

GCMEDICA is a professional and 100% purely manufacturer, who with the whole manufacturing process completed within our plants, from design, manufacture to sterilization.

GCMEDICA upholds “devoted to improve the lives of the people” as its spirit of enterprise. Our commitment to social engagement and responsibility is an important aspect of our company culture and philosophy.

GCMEDICA keeps the concept of being oriented by science, technology and taking quality as life. Providing the reliable products and credit has always developed its advantage. Making the products as a faith, observing market is guides to lead GCMEDICA a world famous brand.




## CONTENTS

|   |       |
|---|-------|
| ◆ UROLOGICAL .....                      | 1~13  |
| ※ CYSTO/BLADDER IRRIGATION SET .....    | 1~5   |
| • M-Easy Bladder Irrigation Set .....   | 1     |
| • B-Cylind Bladder Irrigation Set ..... | 1     |
| • S-Tur Bladder Irrigation Set .....    | 2     |
| • S-Uni Bladder Irrigation Set .....    | 2     |
| • B-Uro Bladder Irrigation Set .....    | 3     |
| • Premi Bladder Irrigation Set .....    | 3     |
| • J-Pump Bladder Irrigation Set .....   | 4     |
| • J-Tur Bladder Irrigation Set .....    | 4     |
| • H-Pump Bladder Irrigation Set .....   | 5     |
| • Sup-Flow Bladder Irrigation Set ..... | 5     |
| ※ ARTHROSCOPY IRRIGATION SET .....      | 6     |
| ※ DRAINAGE BAG WITH URINE METER .....   | 7     |
| ※ LOW BED URINE BAG .....               | 7     |
| ※ URINARY DRAINAGE BAG.....             | 8     |
| ※ ENEMA KITS .....                      | 9     |
| ※ SITZ BATH KITS .....                  | 9     |
| ※ URINARY DRAINAGE LEG BAG .....        | 10    |
| ※ CLICK SEAL SPECIMEN CONTAINER .....   | 11    |
| ※ SUCTION SPECIMEN DRAIN BAG .....      | 11    |
| ※ NELATON CATHETER .....                | 12    |
| ※ SPIGOT CATHETER AND ADAPTOR .....     | 12    |
| ※ SILICON MALE CATHETER .....           | 13    |
| ※ RECTAL BALLOON CATHETER .....         | 13    |
| ◆ GYNAECOLOGY .....                     | 14~16 |
| ※ UTERINE ASPIRATION TUBING .....       | 14    |
| ※ UTERINE ASPIRATION CURETTE .....      | 14    |
| ※ VAGINAL SPECULUM .....                | 15~16 |
| ※ AMNIO HOOK .....                      | 16    |
| ※ IUD REMOVAL HOOK .....                | 16    |

## CYSTO / BLADDER IRRIGATION SET



### M-Easy Bladder Irrigation Set



| Product code         | 081101  | 081102   | 081104  |
|----------------------|---|--|---|
| Non vented spike     | 1   | 2  | 4   |
| Pillar drip chamber  | 1   | 1  | 1   |
| Click clamp          |   | 2  | 4   |
| Y-connector          |   | 1  | 3   |
| Roller regulator     | 1   | 1  | 1   |
| Silicon tube adapter | 1   | 1  | 1   |
| Luer lock connector  |   |  |   |
| Drawing              |  |  |  |





### B-Cylind Bladder Irrigation Set




| Product code         | 081201   | 081202  |
|----------------------|--|---|
| Non vented spike     | 1  | 2   |
| Pillar drip chamber  | 1  | 1   |
| Click clamp          | 2  |   |
| Y-connector          |  |   |
| Roller regulator     | 1  | 1   |
| Silicon tube adapter |  | 1   |
| Drawing              |  |  |

## S-Tur Bladder Irrigation Set




| Product code         | 081501  | 081511  | 081502  | 081512  |
|----------------------|---|---|---|---|
| Non vented spike     | 1   | 1   | 2   | 2   |
| Drip chamber         | 1   | 1   | 1   | 1   |
| Click clamp          |   |   | 2   | 2   |
| Y-connector          |   | 1   |   | 1   |
| Roller regulator     | 1   | 1   | 1   | 1   |
| Silicon tube adapter | 1   | 1   | 1   | 1   |
| Lure lock connector  |   | 1   |   | 1   |
| Drawing              |  |  |  |  |

## S-Uni Bladder Irrigation Set

| Product code         | 081521  |
|----------------------|---|
| Non vented spike     | 1   |
| Drip chamber         | 1   |
| Roller regulator     | 1   |
| Silicon tube adapter | 1   |
| Drawing              |  |




## B-Uro Bladder Irrigation Set



| Product code         | 081401 | 081402 | 081404 |
|----------------------|--------|--------|--------|
| Non vented spike     | 1      | 2      | 4      |
| Drip chamber         | 1      | 1      | 1      |
| Click clamp          |        | 2      | 4      |
| Y-connector          |        | 1      | 3      |
| Roller regulator     | 1      | 1      | 1      |
| Silicon tube adapter | 1      | 1      | 1      |

Drawing



## Premi Bladder Irrigation Set

| Product code         | 082302 | 082304 |
|----------------------|--------|--------|
| Non vented spike     | 2      | 4      |
| Jumbo drip chamber   | 1      | 1      |
| Click clamp          |        |        |
| Y-connector          |        | 3      |
| Roller regulator     |        | 1      |
| Silicon tube adapter | 1      | 1      |

Drawing




## J-Pump Bladder Irrigation Set



| Product code         | 081702 | 081722 |
|----------------------|--------|--------|
| Non vented spike     | 2      | 2      |
| Drip chamber         | 1      | 2      |
| Roller regulator     | 1      | 1      |
| Silicon tube adapter | 1      | 1      |
| Luer lock connector  |        | 1      |

Drawing



## J-Tur Bladder Irrigation Set





| Product code         | 081712 | 081732 |
|----------------------|--------|--------|
| Non vented spike     | 2      | 2      |
| Roller regulator     | 1      | 1      |
| Silicon tube adapter | 1      | 1      |
| Luer lock            |        | 1      |

Drawing



## H-Pump Bladder Irrigation Set




| Product code         | 081801   | 081802   |
|----------------------|--|--|
| Non vented spike     | 1  | 2  |
| Drip chamber         | 1  | 1  |
| Roller regulator     | 1  | 1  |
| Silicon tube adapter | 1  | 1  |
| Luer lock connector  |  | 1  |
| Drawing              |  |  |

## Sup-Flow Bladder Irrigation Set

### FEATURE:





·Burette drip chamber and extra wider tube I.D. for large flow irrigation.



| Product code         | 084404   |
|----------------------|--|
| Non vented spike     | 4  |
| Drip chamber         | 1  |
| Roller regulator     | 1  |
| Silicon tube adapter | 1  |
| Drawing              |  |

## ARTHROSCOPY IRRIGATION SET



| Product code         | 080601  | 080602  | 080611  | 080604  |
|----------------------|---|---|---|---|
| Non vented spike     | 2   | 2   | 2   | 4   |
| Drip chamber         | 1   | 2   | 1   | 1   |
| Click clamp          | 3   | 3   | 3   | 5   |
| Y-connector          | 1   | 1   | 1   | 3   |
| Roller clamp         |   |   | 1   |   |
| Silicon tube adapter | 1   | 1   | 1   | 1   |
| Drawing              |  |  |  |  |

## DRAINAGE BAG WITH URINE METER

Drainage Bag with Urine Meter is a closed system used to accurately measure urine output and collect the urine. The device features a rigid, clear and transparent measured volume chamber, and a 2600mL soft drainage bag.

Product code 051E01

### FEATURES:

- 2600ml drainage bag with 500ml meter
- Anti-reflux valve
- Access port
- Double wrapped
- With convenient hanger
- 150cm tubing
- Latex free
- Single patient use



## LOW BED URINE BAG

Low-Bed Drainage Bag features to facilitate proper drainage and anti-reflux valve for back-flow prevention. The exclusive slide-tap drainage port allows for minimal splashing during drainage. The drainage Bag has a reversible design for versatile placement on either side of the bed. This drain bag hangs only 6-inches (18cm) to allow low placement on an adjustable bed without touching the floor.



Product code 051M01

### FEATURES:

- With cross valve. With anti-reflux valve
- With reinforced double hanger & rope hanger

## URINARY DRAINAGE BAG

- 2000ml,with T-Tap outlet valve,and soft tube pocket
- Needless sample port with cap
- Two parts anti-reflux drip chamber,air-vent on the bag
- With reinforced double hanger & rope hanger
- With bed sheet clamp
- The tubing is  $\varnothing$  10.0×1100mm
- Latex-free,packed sterile

Product code 051B01

### FEATURES:

- Diamond shaped bag
- With cross valve. With anti-reflux valve
- With reinforced single hanger



Product code 051B02

### FEATURES:

- Pear-shaped bag
- With reinforced single hanger & rope hanger

Product code 051D01

### FEATURES:

- Triangle shaped bag
- With single reinforced hanger



## ENEMA KITS

### FEATURES:

- 1500ml bag or bucket with 60 inches tubing, with pre-lubricated tip
- Soap packet
- Water-proof drape: 420x420mm



| Product code | Description           |
|--------------|-----------------------|
| 056A15       | Enema bag only        |
| 056B15       | Enema kit with bucket |

## SITZ BATH KITS

### FEATURES:

- Ideal for the treatment of hemorrhoids or other perineal conditions.
- Provides wide, contoured edges for maximum comfort and support.
- Fits most standard toilets and commodes.
- Water bag can be positioned on toilet seat or on top of toilet tank.
- Easy to clean durable plastic basin.
- 2000ml bag with 60 inches tubing.
- Fits standard toilets and commodes.
- Portable and can administrate soothing fluids to perinea area.



| Product code | Description                |
|--------------|----------------------------|
| 054A01       | Sitz bag only              |
| 059A01       | Sitz bag with a bath basin |

## URINARY DRAINAGE LEG BAG

| Product code | Spec                       |
|--------------|----------------------------|
| 055C06       | 600ml                      |
| 055C16       | 600ml, with extension tube |
| 055C09       | 900ml                      |
| 055C19       | 900ml,with extension tube  |

### FEATURES:

- Square shape,with anti-reflux valve
- Flip style valve drainage outlet
- Two velcro straps
- With twist valve
- With extension tubing



| Product code | Spec  |
|--------------|-------|
| 055D04       | 450ml |
| 055D07       | 750ml |
| 055D09       | 900ml |

### FEATURES:

- Square shape,with anti-reflux valve
- Twist style valve drainage outlet
- With double high elastic belts

| Product code | Spec   |
|--------------|--------|
| 055B05       | 500ml  |
| 055B07       | 700ml  |
| 055B10       | 1000ml |

Urinary Leg Bag is durably constructed for odor containment and leak resistance. The anti-reflux chamber reduces gross reflux.

### FEATURES:

- Elliptic shape , with anti-reflux valve
- Equipped with handle/hanger
- Design allows for even filling and draining
- Includes two flexible straps
- Single use
- Latex free



## CLICK SEAL SPECIMEN CONTAINER

### FEATURES:

- With a sound of click when cap is sealed tight, decreasing the risk of contamination.
- Special cap design ensures convenience and easy job for nurses.
- With graduation for fast volume reading.
- Available with 90ml and 120ml capacity.



| Product code | Description |
|--------------|-------------|
| 061A03       | 90ml        |
| 061A04       | 120ml       |

## SUCTION SPECIMEN DRAIN BAG


### FEATURES:

- Open mesh filter allow for fast drainage of fluids
- Durable, latex drape has flexible hose
- Velcro securement eliminates the need to use tape



| Product code  | 05701            | 05702            |
|---------------|------------------|------------------|
| Filter        | one-piece filter | removable filter |
| Flexible hose | with             | with             |

Drawing



## NELATON CATHETER



| Product code | Description   |
|--------------|---|
| 05801F       | 26cm length tube with smooth two eyes, funnel connector       |
| 05802F       | 26cm length tube with two eyes, funnel connector              |
| 05801M       | 40cm length tube with smooth two eyes, funnel connector       |
| 05802M       | 40cm length tube with two eyes, funnel connector              |
| 058A         | 40cm length tube with one eye, closed tip, connector with cap |

## SPIGOT CATHETER AND ADAPTOR

Catheter Spigot is used to provide stop of flow for catheters hygienically during nursing procedure. It is non-invasive utilized to seal the catheter for a short time to allow urine collect in bladder.

### FEATURES:

- Graduated design makes it easily accommodate tubing of various internal diameters.
- Ribbed design is available.
- Two sizes are available for different tube diameter.



| Product code | Description            |
|--------------|------------------------|
| 1971         | 0.D.10mm               |
| 1972         | 0.D.13mm               |
| 101301       | Ribbed type            |
| 1939A01      | Foley catheter adaptor |




## SILICON MALE CATHETER

### FEATURES:

Silicon male catheter is a self-adhering catheter that is made from odor-free silicone. It is latex-free to reduce skin irritation. Transparent and clear, allowing visual inspections without the need to remove the catheter. Built in adhesive makes applications easier and allows longer wear time. A kink-proof funnel is designed to facilitate uninterrupted urine flow.

- Self-Adhering Bond.
- Kink Proof Funnel for Uninterrupted Flow.
- Latex-Free to Avoid Irritation.
- Constructed with 100% Silicone.
- Transparent for Easy Visual Skin Inspections.
- Hair protector and size indicator is available for convenient application.




| Product code | Description                                       |
|--------------|---|
| 063A25       | I.D. 25mm   |
| 063A29       | I.D. 29mm   |
| 063A32       | I.D. 32mm   |
| 063A36       | I.D. 36mm   |
| 063A41       | I.D. 41mm   |
| 063B25       | I.D. 25mm, With hair protector and size indicator |
| 063B29       | I.D. 29mm, With hair protector and size indicator |
| 063B32       | I.D. 32mm, With hair protector and size indicator |
| 063B36       | I.D. 36mm, With hair protector and size indicator |
| 063B41       | I.D. 41mm, With hair protector and size indicator |

## RECTAL BALLOON CATHETER

### FEATURES:

- Used for measuring rectal pressure during urodynamic studies.
- Silicone construction avoids latex sensitivity concerns.
- Made of PVC with silicone balloon, single lumen
- Universal adapter with sealing cap
- Sterile, single use




|              |         |
|--------------|---------|
| Product code | 6035C01 |
|--------------|---------|

## UTERINE ASPIRATION TUBING

### FEATURES:

- Flexible cuff for convenient operations.
- Special wide bore tubing, thick walled, high vacuum PVC tubing-proof seals.
- Available with flexible cuff, providing extra flexibility at crucial point.
- Easily connects to Curettes.




| Product code | Description  |
|--------------|--|
| 7050A01      | Rotary handle only                                     |
| 7051A01      | I.D.12mm bore with rotary handle                       |
| 7051B01      | I.D.12mm with rotary handle and adapter                |
| 7051C01      | I.D.12mm with rotary handle, adapter and flexible cuff |

## UTERINE ASPIRATION CURETTE

### FEATURES:

- Rigid, clear, plastic curettes for uterine aspiration
- Molded markings for size identification
- Available with straight one and curved one
- Available sizes: 6mm, 7mm, 8mm, 9mm, 10mm, 11mm, 12mm, 14mm, 16mm



| Product code | Description |
|--------------|-------------|
| 7049A        | Straight    |
| 7049B        | Curved      |

## VAGINAL SPECULUM

### FEATURES:


- Crystal clear blade for accurate gynecologic checking
- Unbreakable design for patient's safety.
- Smooth and rounded end for patient safety and comfort.
- Transparent finish yields excellent visibility, sampling and examination accuracy.
- With middle screw type design.



| Product code | Description         |
|--------------|---------------------|
| 7048L01      | Size L              |
| 7048L02      | Size L, unbreakable |
| 7048M01      | Size M              |
| 7048M02      | Size M, unbreakable |
| 7048S01      | Size S              |
| 7048S02      | Size S, unbreakable |

### FEATURES:

- Made from non-toxic material.
- Available with common material and unbreakable.
- Middle screw type.



| Product code | Description          |
|--------------|----------------------|
| 7047L01      | Size L, french type  |
| 7047M01      | Size M, french type  |
| 7047S01      | Size S, french type  |
| 7047XS01     | Size XS, french type |

### FEATURES:

- Made from non-toxic material.
- Available with common material and unbreakable.
- Middle screw type.



| Product code | Description             |
|--------------|-------------------------|
| 7045L01      | Size L, side screw type |
| 7045M01      | Size M, side screw type |
| 7045S01      | Size S, side screw type |

## AMNIO HOOK

### FEATURES:

- Used during vaginal deliveries to puncture the amniotic sac.
- Disposable to prevent infection.



| Product code | Description |
|--------------|-------------|
| 7053A01      |             |

## IUD REMOVAL HOOK

### FEATURES:

- For easy and effective retrieval and removal of an IUD.
- Disposable to prevent infection.



| Product code | Description |
|--------------|-------------|
| 7057A01      |             |